



“Radiant Beauty in Your Hand”

Baby's skin is in the tender state with plenty of moisture, vigorous cell activity and unharmed by the outside environment. As we age, our cell activity decreases, gravity and the environment takes it toll on our skin, causing the repairing ability of the coenzymes and collagen to slow down and even stop, thus resulting in wrinkles.

The **No Lines™ Galvanic Wrinkle Eraser** uses state-of-the-art technology originally developed in Silicon Valley. It uses micro current and ion activity to reverse the visible signs of aging by stimulating collagen in the skin, thereby gradually smoothing the skin and improving the appearance of fine lines, wrinkles, black circles, eye puffiness and crow's feet.

The unique ergonomic shape of the **No Lines™ Galvanic Wrinkle Eraser** is specifically designed to target the delicate areas around the eyes to regain and sustain a more youthful, radiant appearance for life.

The 3 Types of Wrinkles

There are 3 different categories of wrinkles: deep wrinkles, moderate wrinkles and slight wrinkles, based on wrinkles formed at different ages.

Deep wrinkle: formed at about 45-years-old and up.

Moderate wrinkle: formed between 30 to 45 years old

Slight wrinkle: formed below 30-years old.

Scientifically designed to target all 3 wrinkle types (all age levels).

The **No Lines™ Galvanic Wrinkle Eraser** features 3 levels of intensity.

High bionic wave: for smoothing deep wrinkles

Moderate bionic wave: for smoothing moderate wrinkles

Low bionic wave: for smoothing slight wrinkles

Choose the intensity level that corresponds to your needs.

You can gradually decrease the intensity level as the appearance of your wrinkles decrease.



Product Features:

* Power button: On, Off.

* Indicator light of high wave, medium wave, low wave energy
(Adjust to high, medium and low wave energy)

Benefits:

- Minimizes the appearance of fine lines and wrinkles.
- Minimize the appearance of dark circles.
- Slows down the visible effects of aging.
- Designed to target specific areas of concern.
- Immediate and cumulative results.

Operation of the No Lines™ Galvanic Wrinkle Eraser:

1. Press "Power" button for 3-4 seconds to turn the power on.
2. Choose Low, Medium or High level according to your needs by pressing the button.
3. Press "Power" button for 3-4 seconds to turn the power off.
4. Uses one 23A, 12 volt battery (included).

Directions:

- 1) Cleanse with Essona Organics Spring Dew or Seatopia Cleanser. Pat dry.
- 2) Spray Essona SeaMist Tonique on cotton ball and dab on face.
- 3) Apply Essona Organics NuFusion Treatment Serum around the eyes and other target areas.
- 4) Turn on the machine by pressing "Power" switch for 3-4 seconds
- 5) Choose desired intensity level of wave energy by pressing "Power" switch, High, Medium or Low
- 6) Hold the unit comfortably as you would a pen with fingers on shiny sensor area.
- 7) Contact the tip to the wrinkles, indicator flashes and buzzer beeps. The machine is working properly. Move it slowly from the outer to inner face areas, around the eye area. Focus on one area (wrinkle) for 5-10 seconds. Treat four times for every individual wrinkle to achieve best results.
- 8) Turn the power off by pressing "Power" switch for 3-4 seconds
- 9) Apply Essona NoLines Organic Eye Serum gently with the fingertips.

Notes:

- * Use 2-3 times a week. Can be used daily for faster results.
- * Wait at least 6 hours between uses.
- * Keep skin around eyes clean and moisturized.
- * We suggest taking pictures of your eyes in close-up at the same time every 1-2 weeks to enjoy the rejuvenating results.

Cautions:

- * Do not disassemble the unit. No user serviceable parts inside.
 - * Do not drop.
 - * Avoid overexposure to the sun.
 - * Do not use while bathing or showering.
- Do not use if you have disease of the eyes.

**The No Lines™ Galvanic Wrinkle Eraser
is designed to be used with Essona Organics Skin Care products.**



Inner Health, Outer Beauty

www.essona.com